

Using Microscopes


Total Magnification = magnification of OBJECTIVE lens x magnification of EYEPIECE lens

Rules for microscope drawing:

1. Use a sharp pencil
2. Don't colour/shade in any parts of the drawing
3. Give your drawing a heading
4. Draw the specimen in a circle
5. Calculate the total magnification and write it on your drawing


Practice drawing microscope diagrams:

Drawing of _____




Total Magnification = _____ x _____
= _____

Drawing of _____



Total Magnification = _____ x _____
= _____

Drawing of _____



Total Magnification = _____ x _____
= _____